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Version 1.00

Oscilla[®]

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

Classification:

Group 1, class B EN 60601-1-2:2002

Type of protection against electric shock:	Class I equipment
Degree of protection against electric shock:	Type B applied part
Degree of protection against liquid penetration:	IPO, ordinary equipment
Degree of safety of application in the presence of flammable anaesthetics:	N/A
Mode of operation:	Continuous operation

Medical CE- mark:

Inmedico A/S is approved for medical CE marking, by DGM.
Identification number 0543

Transducers:

TDH-39 air

Power Supply:

12-18 volt AC or DC, 15 watt max.

Adaptor:

Primary: 230 VAC, 50/60 Hz
Secondary: 18 VAC, 20 VA, $\pm 10\%$

Environmental Conditions for Operation

Ambient Temperature: +15 to +35 degree Celsius.
Relative Humidity: 30 % to 90 %
Surrounding pressure: 80 kPa to 120 kPa

Environmental Conditions for Storage

Ambient Temperature: -10 to +50 degree Celsius
Relative Humidity: 95% or less (non-condensing)
Surroundings pressure: 50 kPa to 120 kPa

Physical Attributes

Dimensions: 295 (W) x 180 (D) x 55 (H) mm (front : 25 mm)
Weight: ~ 0,9 kg including accessories
Casing material: Aluminium

Calibration:

Data is saved in eeprom memory, EN 60604-1

Warm-up time:

< 3 seconds.

Included parts:

1 audiometer, 1 set of Peltor noise reducing headphones mounted with telephones,
1 response button, 1 mains adapter, 1 user manual and calibration certificate.

Accessories: carrying bag, serial cable, AudioConsole software

Measurement specifications

Measurement Method Manual and automatic threshold test **with** storage of data.

Measurement Range
Maximum Intensities:

Frequency	Air
125 Hz	70 dB
250 Hz	90 dB
500 Hz	110 dB
750 Hz	110 dB
1000 Hz	110 dB
1500 Hz	110 dB
2000 Hz	110 dB
3000 Hz	110 dB
4000 Hz	110 dB
6000 Hz	100 dB
8000 Hz	90 dB

Measurement Accuracy

Tolerance frequency: $\pm 1\%$.
In/output impedance: Left / Right channel $< 1\text{ ohm}$, Patient response $\sim 10\text{ K ohm}$
Distortion speaker: TDH-39: $< 1\%$, 4. and higher harmonic + sub harmonic $< 0,3\%$.

Single signal channel:
Air conductor: $\pm 4\text{ dB}$ (125 Hz – 4 kHz)
 $\pm 5\text{ dB}$ (6 kHz – 8 kHz)

Noise reduction

Peltor H7A (Optime II): <http://www.peltor.se>

Frequency	125	250	500	1000	2000	4000	8000	H	M	L	SNR
Mean att.	14.5	20.3	32.6	39.1	35.1	34.7	39.8				
Stand. dev.	1.8	2.6	2.4	2.5	2.3	2.7	2.5	34 DB	28 DB	20 DB	31 DB
APV-value	12.7	17.7	30.2	36.6	32.8	32.0	37.3				

SNR (Single Number Rating or simplified Noise Level Reduction) indicates the attenuation level in dB in an average industrial noise environment.

(H), (M) and (L) values in dB indicate the attenuation in high frequency (H), medium frequency (M) and in low frequency (L) noise.

6 Please remark

6.1 Safety precaution

Medical electrical equipment need special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the accompanying documents.

Portable and mobile RF communications equipment can affect medical electrical equipment.

The audiometer should be placed in such a way that the patient can not see or hear when it is operated. If a sound booth is not used, place the audiometer at least 1 meter from the patient. Operate the interrupter buttons with a light soft touch to prevent the patient from hearing any sounds from the buttons.

Interruptions of the mains adapter must be carried out on the outlet.

If the audiometer / personal computer are exposed to a powerful static discharge, this can cause interruption of the signal without interruption of the device.

Should this rare state occur? Turn the audiometer off and then on. Under normal conditions the device should function again.

Neither the patient nor user will be exposed to any danger; likewise the calibrations settings or other functions would be affected.

If a failure should occur, meaning the device is missing a function or has a fault, please contact the seller of the device or the manufacture.

Repair of the device must be accomplished by Inmedico A/S or special trained persons with a given permission to repair the product.

PC, printers or other equipment there are used together with the audiometer, must meet the specified demands there are given in standard EN60601-1 § 14.1a

Always use accessories, such as, attachments or spare parts, provided for exclusive use of this instrument. Using other than them may cause the instrument to fail or fulfill its function completely.

6.2 Packing guidance

By shipment, the device must be wrapped into bubble plastic and put into an ordinary export case.

None of the audiometer parts must be able to touch the sides in the box.

The audiometer withstand ordinary air freight and similar.

6.3 Calibration

It's recommended to get the device Calibrated every 2nd year, by qualified personal.

6.4 Cleaning / disinfection

Audiometer parts that are in contact with the patient include patient response and headset.

Both can be cleaned with a cloth wrung in lukewarm water and a little soap. Disconnect both parts from the audiometer before cleaning and make sure they are dry before connecting them again.

6.5 Product Warranty

2 years – device only

6.6 Product Life

20 years - estimated

7 EMC

Guidance and manufacture's declaration – electromagnetic emission

The SM920-P is intended for use in the electromagnetic environment specified below. The customer or the user of the SM920-P should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The SM920-P uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The SM920-P is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Complies	
fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacture's declaration – electromagnetic immunity

The SM920-P is intended for use in the electromagnetic environment specified below. The customer or the user of the SM920-P should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-5	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycle 70 % U_T (30 % dip in U_T) for 25 cycle <5 % U_T (>95 % dip in U_T) for 5 sec	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycle 70 % U_T (30 % dip in U_T) for 25 cycle <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the SM920-P requires continued operation during power mains interruptions, it is recommended that the SM920-P be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacture's declaration – electromagnetic immunity

The SM920-P is intended for use in the electromagnetic environment specified below. The customer or the user of the SM920-P should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the SM920-P, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter Recommended separation distance $d = 1,17\sqrt{P}$ $d = 1,17\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d = 2,33\sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ MHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

CALIBRATION CERTIFICATE FOR AUDIOMETER

Serial no.: 920P-1406-AUDF

Type: Oscilla SM920-P Screening audiometer

Customer: Audifon Mróz Sp.J - Poland

Firmware: 1.01

Print: 28173-H-20

Phone1 - Transducer: TDH-39

Reference standard: EN60645

Serial. Left: D009184

Serial. Right: D009202

Frequency Hz	RETSPL ISO 389-1 1998	Test dBSPL	Accepted Deviation dB	Soundpressure Tone		Frequency Accepted Deviation%
				Left	Right	
125	45	80	+/- 3	✓	✓	✓ +/- 1
250	25,5	80	+/- 3	✓	✓	✓ +/- 1
500	11,5	80	+/- 3	✓	✓	✓ +/- 1
750	7,5	80	+/- 3	✓	✓	✓ +/- 1
1000	7	80	+/- 3	✓	✓	✓ +/- 1
1500	6,5	80	+/- 3	✓	✓	✓ +/- 1
2000	9	80	+/- 3	✓	✓	✓ +/- 1
3000	10	80	+/- 3	✓	✓	✓ +/- 1
4000	9,5	80	+/- 3	✓	✓	✓ +/- 1
6000	15,5	80	+/- 5	✓	✓	✓ +/- 1
8000	13	80	+/- 5	✓	✓	✓ +/- 1

* : IEC 60318-3 (6cc coupler)

Levelmeter: BK3560 SN:2630546

Artificial Ear: BK4152 SN:994620

Date:

Sign.:

Authorized sign.:

2014-12-30

Theresa T. Amato

Jan B. Reel